



EGFL7 Antibody, Biotin conjugated

Product Code	CSB-PA890692LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHF1
Immunogen	Recombinant Human Epidermal growth factor-like protein 7 protein (73-198AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Epidermal growth factor-like protein 7 (EGF-like protein 7) (Multiple epidermal growth factor-like domains protein 7) (Multiple EGF-like domains protein 7) (NOTCH4-like protein) (Vascular endothelial statin) (VE-statin) (Zneu1), EGFL7, MEGF7
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	EGFL7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.