



FBXO10 Antibody, HRP conjugated

Product Code	CSB-PA890915LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UK96
Immunogen	Recombinant Human F-box only protein 10 protein (248-400AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-box only protein 10, FBXO10, FBX10 PRMT11
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FBXO10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.