




SLC45A2 Antibody

Product Code	CSB-PA890941LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UMX9
Immunogen	Recombinant Human Membrane-associated transporter protein (1-46AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Membrane-associated transporter protein (Melanoma antigen AIM1) (Protein AIM-1) (Solute carrier family 45 member 2), SLC45A2, AIM1 MATP
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	SLC45A2

Image



Western Blot

Positive WB detected in: MCF-7 whole cell lysate, HeLa whole cell lysate

All lanes: SLC45A2 antibody at 5.3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 59, 34, 52 kDa

Observed band size: 59 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.