



SLC45A2 Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA890941LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9UMX9 |
| Immunogen | Recombinant Human Membrane-associated transporter protein (1-46AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Membrane-associated transporter protein (Melanoma antigen AIM1) (Protein AIM-1) (Solute carrier family 45 member 2), SLC45A2, AIM1 MATP |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Tags & Cell Markers |
| Target Names | SLC45A2 |