



RAB21 Antibody, Biotin conjugated

Product Code	CSB-PA891733LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UL25
Immunogen	Recombinant Human Ras-related protein Rab-21 protein (141-222AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-related protein Rab-21, RAB21, KIAA0118
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RAB21
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.