



PLAGL2 Antibody, HRP conjugated

Product Code	CSB-PA891980LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UPG8
Immunogen	Recombinant Human Zinc finger protein PLAGL2 protein (324-469AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein PLAGL2 (Pleiomorphic adenoma-like protein 2), PLAGL2, KIAA0198
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PLAGL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.