



ABCC6 Antibody, Biotin conjugated

Product Code	CSB-PA001065LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95255
Immunogen	Recombinant Human Multidrug resistance-associated protein 6 protein (730-931AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Multidrug resistance-associated protein 6 (ATP-binding cassette sub-family C member 6) (Anthracycline resistance-associated protein) (Multi-specific organic anion transporter E) (MOAT-E), ABCC6, ARA MRP6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ABCC6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.