



TYW5 Antibody, Biotin conjugated

Product Code	CSB-PA003860LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A2RUC4
Immunogen	Recombinant Human tRNA wybutosine-synthesizing protein 5 protein (1-152AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	tRNA wybutosine-synthesizing protein 5 (hTYW5) (EC 1.14.11.42) (tRNA(Phe) (7-(3-amino-3-carboxypropyl)wyosine(37)-C(2))-hydroxylase), TYW5, C2orf60
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TYW5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.