



CCT6A Antibody, HRP conjugated

Product Code	CSB-PA004865EB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40227
Immunogen	Recombinant Human T-complex protein 1 subunit zeta protein (80-250AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-complex protein 1 subunit zeta (TCP-1-zeta) (Acute morphine dependence-related protein 2) (CCT-zeta-1) (HTR3) (Tcp20), CCT6A, CCT6 CCTZ
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CCT6A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.