




cyp26a1 Antibody

Product Code	CSB-PA006402LA01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P79739
Immunogen	Recombinant Zebrafish Cytochrome P450 26A1 protein (244-371AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytochrome P450 26A1 (EC 1.14.13.-) (Cytochrome P450RAI) (Retinoic acid 4-hydroxylase) (Retinoic acid-metabolizing cytochrome), cyp26a1, cyp26
Immunogen Species	Zebrafish
Research Area	Others
Target Names	cyp26a1

Image



Western Blot

Positive WB detected in: Zebrafish tissue, Zebrafish eye tissue

All lanes: cyp26a1 at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 57 kDa

Observed band size: 57 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.