



DHCR7 Antibody, Biotin conjugated

Product Code	CSB-PA006844LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBM7
Immunogen	Recombinant Human 7-dehydrocholesterol reductase protein (63-147AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	7-dehydrocholesterol reductase (7-DHC reductase) (EC 1.3.1.21) (Putative sterol reductase SR-2) (Sterol Delta(7)-reductase), DHCR7, D7SR
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	DHCR7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.