



FOXA3 Antibody, Biotin conjugated

Product Code	CSB-PA008796LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55318
Immunogen	Recombinant Human Hepatocyte nuclear factor 3-gamma protein (220-336AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hepatocyte nuclear factor 3-gamma (HNF-3-gamma) (HNF-3G) (Fork head-related protein FKH H3) (Forkhead box protein A3) (Transcription factor 3G) (TCF-3G), FOXA3, HNF3G TCF3G
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	FOXA3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.