



GLRB Antibody, FITC conjugated

Product Code	CSB-PA009520LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48167
Immunogen	Recombinant Human Glycine receptor subunit beta protein (355-472AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glycine receptor subunit beta (Glycine receptor 58 kDa subunit), GLRB
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GLRB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.