



GUCY2D Antibody, HRP conjugated

Product Code	CSB-PA010054LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02846
Immunogen	Recombinant Human Retinal guanylyl cyclase 1 protein (362-462AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Retinal guanylyl cyclase 1 (RETGC-1) (EC 4.6.1.2) (Guanylate cyclase 2D, retinal) (Rod outer segment membrane guanylate cyclase) (ROS-GC), GUCY2D, CORD6 GUC1A4 GUC2D RETGC RETGC1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GUCY2D
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.