



HBG2 Antibody, Biotin conjugated

Product Code	CSB-PA010156LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P69892
Immunogen	Recombinant Human Hemoglobin subunit gamma-2 protein (2-147AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hemoglobin subunit gamma-2 (Gamma-2-globin) (Hb F Ggamma) (Hemoglobin gamma-2 chain) (Hemoglobin gamma-G chain), HBG2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	HBG2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.