



HOXB3 Antibody, Biotin conjugated

Product Code	CSB-PA010663LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14651
Immunogen	Recombinant Human Homeobox protein Hox-B3 protein (315-423AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox protein Hox-B3 (Homeobox protein Hox-2.7) (Homeobox protein Hox-2G), HOXB3, HOX2G
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HOXB3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.