



IGHA2 Antibody, Biotin conjugated

Product Code	CSB-PA011112LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01877
Immunogen	Recombinant Human Immunoglobulin heavy constant alpha 2 protein (36-111AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Immunoglobulin heavy constant alpha 2 (Ig alpha-2 chain C region) (Ig alpha-2 chain C region BUT) (Ig alpha-2 chain C region LAN), IGH2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IGHA2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.