



IRX6 Antibody, HRP conjugated

Product Code	CSB-PA011838LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78412
Immunogen	Recombinant Human Iroquois-class homeodomain protein IRX-6 protein (289-389AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Iroquois-class homeodomain protein IRX-6 (Homeodomain protein IRXB3) (Iroquois homeobox protein 6), IRX6, IRX7 IRXB3
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	IRX6