



IRX6 Antibody, Biotin conjugated

Product Code	CSB-PA011838LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78412
Immunogen	Recombinant Human Iroquois-class homeodomain protein IRX-6 protein (289-389AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Iroquois-class homeodomain protein IRX-6 (Homeodomain protein IRXB3) (Iroquois homeobox protein 6), IRX6, IRX7 IRXB3
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	IRX6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.