



KIF3C Antibody, HRP conjugated

Product Code	CSB-PA012337LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14782
Immunogen	Recombinant Human Kinesin-like protein KIF3C protein (158-267AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kinesin-like protein KIF3C, KIF3C
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KIF3C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.