



KRIT1 Antibody, Biotin conjugated

Product Code	CSB-PA012501LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00522
Immunogen	Recombinant Human Krev interaction trapped protein 1 protein (471-585AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Krev interaction trapped protein 1 (Krev interaction trapped 1) (Cerebral cavernous malformations 1 protein), KRIT1, CCM1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KRIT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.