



LILRB3 Antibody, Biotin conjugated

Product Code	CSB-PA012941LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75022
Immunogen	Recombinant Human Leukocyte immunoglobulin-like receptor subfamily B member 3 protein (470-624AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leukocyte immunoglobulin-like receptor subfamily B member 3 (LIR-3) (Leukocyte immunoglobulin-like receptor 3) (CD85 antigen-like family member A) (Immunoglobulin-like transcript 5) (ILT-5) (Monocyte inhibitory receptor HL9) (CD antigen CD85a), LILRB3, ILT5 LIR3
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	LILRB3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.