



LPO Antibody, HRP conjugated

Product Code	CSB-PA013066LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22079
Immunogen	Recombinant Human Lactoperoxidase protein (81-190AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lactoperoxidase (LPO) (EC 1.11.1.7) (Salivary peroxidase) (SPO), LPO, SAPX
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	LPO
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.