



LPO Antibody, FITC conjugated

Product Code	CSB-PA013066LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22079
Immunogen	Recombinant Human Lactoperoxidase protein (81-190AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lactoperoxidase (LPO) (EC 1.11.1.7) (Salivary peroxidase) (SPO), LPO, SAPX
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	LPO