



LRIG2 Antibody, Biotin conjugated

| | |
|---------------------|--|
| Product Code | CSB-PA013085LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O94898 |
| Immunogen | Recombinant Human Leucine-rich repeats and immunoglobulin-like domains protein 2 protein (833-932AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Leucine-rich repeats and immunoglobulin-like domains protein 2 (LIG-2), LRIG2, KIAA0806 LIG2 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cancer |
| Target Names | LRIG2 |