



MS4A2 Antibody, Biotin conjugated

Product Code	CSB-PA015013LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01362
Immunogen	Recombinant Human High affinity immunoglobulin epsilon receptor subunit beta protein (1-59AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	High affinity immunoglobulin epsilon receptor subunit beta (FcεRI) (Fc epsilon receptor I beta-chain) (IgE Fc receptor subunit beta) (Membrane-spanning 4-domains subfamily A member 2), MS4A2, APY FCER1B IGER
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MS4A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.