



MT-CO2 Antibody, FITC conjugated

Product Code	CSB-PA015073LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00403
Immunogen	Recombinant Human Cytochrome c oxidase subunit 2 protein (83-227AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytochrome c oxidase subunit 2 (Cytochrome c oxidase polypeptide II), MT-CO2, COII COX2 COXII MTCO2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MT-CO2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.