



MTSS1 Antibody, Biotin conjugated

Product Code	CSB-PA015201LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43312
Immunogen	Recombinant Human Metastasis suppressor protein 1 protein (616-730AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Metastasis suppressor protein 1 (Metastasis suppressor YGL-1) (Missing in metastasis protein), MTSS1, KIAA0429 MIM
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MTSS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.