



NAP1L5 Antibody, FITC conjugated

Product Code	CSB-PA015445LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96NT1
Immunogen	Recombinant Human Nucleosome assembly protein 1-like 5 protein (1-67AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleosome assembly protein 1-like 5 (Down-regulated in liver malignancy), NAP1L5, DRLM
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NAP1L5