



OPN1MW Antibody, Biotin conjugated

Product Code	CSB-PA016352LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04001
Immunogen	Recombinant Human Medium-wave-sensitive opsin 1 protein (1-52AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Medium-wave-sensitive opsin 1 (Green cone photoreceptor pigment) (Green-sensitive opsin) (GOP), OPN1MW, GCP
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	OPN1MW
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.