




PACS1 Antibody

Product Code	CSB-PA017371LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6VY07
Immunogen	Recombinant Human Phosphofurin acidic cluster sorting protein 1 protein (757-905AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphofurin acidic cluster sorting protein 1 (PACS-1), PACS1, KIAA1175
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	PACS1

Image



Western Blot

Positive WB detected in: Raji whole cell lysate,

HeLa whole cell lysate

All lanes: PACS1 antibody at 3.2µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 105, 110 kDa

Observed band size: 110 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.