



PACS1 Antibody, Biotin conjugated

Product Code	CSB-PA017371LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6VY07
Immunogen	Recombinant Human Phosphofurin acidic cluster sorting protein 1 protein (757-905AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphofurin acidic cluster sorting protein 1 (PACS-1), PACS1, KIAA1175
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	PACS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.