



# PARP3 Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA017467EB01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9Y6F1
<b>Immunogen</b>	Recombinant Human Poly [ADP-ribose] polymerase 3 protein (1-240AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Poly [ADP-ribose] polymerase 3 (PARP-3) (hPARP-3) (EC 2.4.2.30) (ADP-ribosyltransferase diphtheria toxin-like 3) (ARTD3) (IRT1) (NAD(+) ADP-ribosyltransferase 3) (ADPRT-3) (Poly[ADP-ribose] synthase 3) (pADPRT-3), PARP3, ADPRT3 ADPRTL3
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	PARP3
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.