



PHEX Antibody, HRP conjugated

Product Code	CSB-PA017896LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78562
Immunogen	Recombinant Human Phosphate-regulating neutral endopeptidase protein (524-673AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphate-regulating neutral endopeptidase (EC 3.4.24.-) (Metalloendopeptidase homolog PEX) (Vitamin D-resistant hypophosphatemic rickets protein) (X-linked hypophosphatemia protein) (HYP), PHEX, PEX
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	PHEX
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.