



# PHEX Antibody, Biotin conjugated

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA017896LD01HU   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| <b>Uniprot No.</b>         | P78562   |
| <b>Immunogen</b>           | Recombinant Human Phosphate-regulating neutral endopeptidase protein (524-673AA)   |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Human  |
| <b>Tested Applications</b> | ELISA  |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | Biotin   |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4   |
| <b>Purification Method</b> | >95%, Protein G purified   |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Alias</b>               | Phosphate-regulating neutral endopeptidase (EC 3.4.24.-)<br>(Metalloendopeptidase homolog PEX) (Vitamin D-resistant hypophosphatemic rickets protein) (X-linked hypophosphatemia protein) (HYP), PHEX, PEX |
| <b>Immunogen Species</b>   | Homo sapiens (Human)   |
| <b>Research Area</b>       | Stem Cells   |
| <b>Target Names</b>        | PHEX   |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.  |