



PLEKHG5 Antibody, FITC conjugated

Product Code	CSB-PA018174LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94827
Immunogen	Recombinant Human Pleckstrin homology domain-containing family G member 5 protein (821-957AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pleckstrin homology domain-containing family G member 5 (PH domain-containing family G member 5) (Guanine nucleotide exchange factor 720) (GEF720), PLEKHG5, KIAA0720
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PLEKHG5