



PLEKHG5 Antibody, Biotin conjugated

Product Code	CSB-PA018174LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94827
Immunogen	Recombinant Human Pleckstrin homology domain-containing family G member 5 protein (821-957AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pleckstrin homology domain-containing family G member 5 (PH domain-containing family G member 5) (Guanine nucleotide exchange factor 720) (GEF720), PLEKHG5, KIAA0720
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PLEKHG5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.