



PLXNC1 Antibody, HRP conjugated

Product Code	CSB-PA018225LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60486
Immunogen	Recombinant Human Plexin-C1 protein (81-218AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Plexin-C1 (Virus-encoded semaphorin protein receptor) (CD antigen CD232), PLXNC1, VESPR
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	PLXNC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.