



PLXNC1 Antibody, FITC conjugated

Product Code	CSB-PA018225LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60486
Immunogen	Recombinant Human Plexin-C1 protein (81-218AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Plexin-C1 (Virus-encoded semaphorin protein receptor) (CD antigen CD232), PLXNC1, VESPR
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	PLXNC1