



POU1F1 Antibody, HRP conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA018391LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P28069 |
| Immunogen | Recombinant Human Pituitary-specific positive transcription factor 1 protein (34-130AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Pituitary-specific positive transcription factor 1 (PIT-1) (Growth hormone factor 1) (GHF-1), POU1F1, GHF1 PIT1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Epigenetics and Nuclear Signaling |
| Target Names | POU1F1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |