



# POU1F1 Antibody, FITC conjugated

Product Code	CSB-PA018391LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28069
Immunogen	Recombinant Human Pituitary-specific positive transcription factor 1 protein (34-130AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pituitary-specific positive transcription factor 1 (PIT-1) (Growth hormone factor 1) (GHF-1), POU1F1, GHF1 PIT1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	POU1F1