



POU3F4 Antibody, FITC conjugated

Product Code	CSB-PA018399LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49335
Immunogen	Recombinant Human POU domain, class 3, transcription factor 4 protein (22-134AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	POU domain, class 3, transcription factor 4 (Brain-specific homeobox/POU domain protein 4) (Brain-4) (Brn-4) (Octamer-binding protein 9) (Oct-9) (Octamer-binding transcription factor 9) (OTF-9), POU3F4, BRN4 OTF9
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	POU3F4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.