



PRB1 Antibody, Biotin conjugated

Product Code	CSB-PA018629LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04280
Immunogen	Recombinant Human Basic salivary proline-rich protein 1 protein (17-92AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Basic salivary proline-rich protein 1 (Salivary proline-rich protein) [Cleaved into: Proline-rich peptide II-2; Basic peptide IB-6; Peptide P-H], PRB1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PRB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.