



PTPN9 Antibody, Biotin conjugated

Product Code	CSB-PA019045LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43378
Immunogen	Recombinant Human Tyrosine-protein phosphatase non-receptor type 9 protein (221-339AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tyrosine-protein phosphatase non-receptor type 9 (EC 3.1.3.48) (Protein-tyrosine phosphatase MEG2) (PTPase MEG2), PTPN9
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PTPN9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.