



QSOX1 Antibody, HRP conjugated

Product Code	CSB-PA019144LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00391
Immunogen	Recombinant Human Sulfhydryl oxidase 1 protein (101-175AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sulfhydryl oxidase 1 (hQSOX) (EC 1.8.3.2) (Quiescin Q6), QSOX1, QSCN6
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	QSOX1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.