



RAB40AL Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA019203LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P0C0E4 |
| Immunogen | Recombinant Human Ras-related protein Rab-40A-like protein (1-88AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Ras-related protein Rab-40A-like (Ras-like GTPase), RAB40AL, RLGP |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | RAB40AL |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |