



RAB40AL Antibody, Biotin conjugated

Product Code	CSB-PA019203LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C0E4
Immunogen	Recombinant Human Ras-related protein Rab-40A-like protein (1-88AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-related protein Rab-40A-like (Ras-like GTPase), RAB40AL, RLGP
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RAB40AL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.