



SERF1A Antibody

Product Code	CSB-PA021040LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75920
Immunogen	Recombinant Human Small EDRK-rich factor 1 protein (1-74AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA;
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Developmental Biology
Target Names	SERF1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.