



SERF1A Antibody, FITC conjugated

Product Code	CSB-PA021040LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75920
Immunogen	Recombinant Human Small EDRK-rich factor 1 protein (1-74AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Small EDRK-rich factor 1 (Protein 4F5) (h4F5) (SMA modifier 1), SERF1A; SERF1B, FAM2A SERF1 SMAM1; FAM2B SERF1 SMAM1
Immunogen Species	Homo sapiens (Human)
Research Area	Developmental Biology
Target Names	SERF1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.