



SLC26A2 Antibody, Biotin conjugated

Product Code	CSB-PA021525LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50443
Immunogen	Recombinant Human Sulfate transporter protein (614-739AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sulfate transporter (Diastrophic dysplasia protein) (Solute carrier family 26 member 2), SLC26A2, DTD DTDST
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC26A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.